

REPLACEMENT PARTS LIST
STERILGARD® II MODELS SG400ULT & SG600ULT
VERTICAL LAMINAR FLOW CABINET

<i>Part Name</i>	<i>SG400ULT Part Numbers BOM SG400ULT#194</i>	<i>SG600ULT Part Numbers BOM SG600ULT#194</i>
Air Pressure Monitor	34697	34697
Capacitor	11557	11558
Circuit Breaker - 5A	34331	34331
Exhaust Filter	37206 (ULPA)	36262 (HEPA)
Fluorescent Ballast	34160	34160
Fluorescent Lamp	17927	17927
Ground Fault Interrupter with Indicator Light	34921	34921
Lamp Holder, "Butt-On" Type	34465	34465
Magnetic Sensor	34609	34609
Membrane Switch Controller	34699	34699
Motor/Blower Controller	304227	304227
Relay, Solid State	33985	33985
Sash Balance - Left Piston	36475	36477
Sash Balance - Right Piston	36476	36478
Supply Blower	11429	11432
Supply Filter	37173 (ULPA)	28303
Supply Motor	31676	31679
Switch, Touch-Type	34615	34615
Transformer	34326	34326
Ultraviolet Lamp	18023	18025
Ultraviolet Lamp Ballast	34160	34165
Ultraviolet Lamp Socket	18088	20281 & 20283
Window Assembly	1851703	1861692
Window Wiper	1855095	1865095

NOTE: When ordering replacement parts, please furnish the serial number of the unit, as well as the model number.