

REPLACEMENT PARTS LIST

PART NAME	QTY Req'd	PART NUMBER
Blower	2	11410
Capacitor	2	11558
Motor (220 v, ½ HP, 1625 RPM)	2	11629
Filter, exhaust, 19" x 38" x 12"	2	39739
Filter, supply, 32" x 45" x 3"	2	38913
Transformer, 120 to 24VAC	1	19913
Motor speed control, 220 V	2	304302
Fluorescent lamp	4	17927
Ballast, fluorescent	2	34160
Differential pressure monitor	1	
Hydraulic pump kit (cart)	1	38726

NOTE: When ordering replacement parts, please furnish serial number of unit, as well as model number