

**REPLACEMENT PARTS FOR  
THE NCB-C4 & NCB-C6**

<b>Part Name</b>	<b>NCB-C4</b>	<b>NCB-C6</b>
Blower	11340	11429
Blower motor	31676	31676
Buzzer	16634	16634
Speed control	304227	304227
Capacitor	11557	11557
Blower switch	15826	15826
Fluorescent light switch	15825	15825
Lamp holders	18123, 18124	18123, 18124
Fluorescent lamp	17966	17968
Proximity sensor	20258	20258
Fluorescent ballast	18168	18167
Circuit breaker switch	186118	186118
HEPA exhaust filter	12114	12117
HEPA supply filter	21848	21951
Work area HEPA filter	12095	12096
Air petcock	19326	19326
Drain valve	12846	12846
View screen glass	32797	17784
View screen cable assembly (RH)	187040A	188086
View screen cable assembly (LH)	187038A	188087
Silicone glass wiper (set)	187074A	188198
AFM monitor	21901	21901
ALM switch	18364	18364

**Optional**

Ultraviolet lamp	18023	18023
Ultraviolet starter	20297	20297
Ultraviolet ballast	18135	18135
Ultraviolet lamp holders	20287	20287
Ultraviolet-fluorescent light switch	15827	15827
Ultraviolet starter base	18082	18082

**NOTE:** When ordering replacement parts, please furnish model number and serial number of unit.