

REPLACEMENT PARTS LIST
STERILGARD III ADVANCE , MODEL SG403A & SG603A
VERTICAL LAMINAR FLOW CABINET

| | SG403A | SG603A |
|--|---------------|---------------|
| Part Description | Part # | Part # |
| Capacitor | 11557 | 11558 |
| Circuit Breaker (Outlets) | 34331 | 34337 |
| Exhaust Hepa Filter | 40661 | 40662 |
| Fluorescent Ballast | 34160 | 40457 |
| Fluorescent Lamp | 17927 | 17927 |
| Front Panel | 333X604 | 336X604 |
| Ground Fault Interrupter Outlet with Indicator Light | 34921 | 34921 |
| Lamp Holder Base | 34465 | 34465 |
| Proximity Switch (Lower) | 40647 | 40647 |
| Proximity Switch (Higher) | 40671 | 40671 |
| Right Angle Micro Cable | 40672 | 40672 |
| Mag Gauge | 20773 | 20773 |
| Membrane Switch Controller | 40702 | 40702 |
| Motor/Blower Speed Controller | 39658 | 39658 |
| Relay, Solid State | 1855167 | 1855167 |
| Terminal Block Asm. | 40450 | 40450 |
| Circuit Breaker (Resettable) | 40449 | 40449 |
| Sash Balance | 38014 | 38172 / 38173 |
| Side Panel, LH (POBM) | 333X650 | 333X650 |
| Side Panel, RH (POBM) | 333X651 | 333X651 |
| Standard Outlet | 18231 | 18231 |
| Supply Blower | 11429 | 11432 |
| Supply Diffuser | 333P709 | 336P709 |
| Supply Hepa Filter | 40141 | 40142 |
| Supply Motor | 31676 | 31679 |
| Switch, Touch Type | 37919 | 37919 |
| Transformer | 34327 | 34327 |
| Ultraviolet Lamp | 18024 | 18025 |
| Ultraviolet Lamp Ballast | 34160 | 34165 |
| Ultraviolet Lamp Socket | 20281 / 20283 | 20281 / 20283 |
| Viewscreen Assembly | 333A310 | 336A310 |
| Viewscreen Frame Assembly | 333A311 | 336A311 |