

Replacement Parts

<i>MODEL</i>	<i>SG403EN</i>	<i>SG603EN</i>
<i>Part Name</i>	<i>Part No. SG403EN</i>	<i>Part No. SG603EN</i>
Capacitor	35617	36517
Circuit Breaker - 5A	34331	34331
Exhaust HEPA Filter	37868	38585
Fluorescent Ballast	34524	34524
Fluorescent Lamp	17927	17927
Front Panel, Dress Panel 10 Degree	333165	336165
Schuko Outlets	39740	39740
Lamp Holder, "Butt-On" Type	34465	34465
Limit Switch, SPDT	16651	16651
Modus, Pressure Monitor	38997	38997
Membrane Switch Controller	34699B	34699B
Motor/Blower Controller	304227	304227
Relay, Solid State	1855167	1855167
Sash Balance	38014	
Sash Balance R.H.		38172
Sash Balance L.H.		38173
Side Panel, LH	333063	333063
Side Panel, RH	333064	333064
Hour Meter	36268	36268
Supply Blower DD10-4AT	11416	
Supply Blower DD10-8AT		11438
Supply Diffuser	333446	336097
Supply HEPA Filter	37865	37866
Supply Motor 1/2 HP 1625RPM	11629	11629
RFI Filter	39630	39630
Switch, Touch-Type	37919	37919
Label, Touchpad, overlay	37998	37998
Transformer	34327	34327
Ultraviolet Lamp	18024	18025
Ultraviolet Lamp Ballast	35421	38241
Ultraviolet Lamp Socket (plunger type)	20281	20281
Ultraviolet Lamp Socket (fixed type)	20283	20283
Viewscreen Assembly	333A300	336A300
Viewscreen Frame Assembly	333081	336081
DC Power Supply	39670	39670
Timer	39820	39820
Limit Switch, DPDT	39252	39252