

<i>MODEL</i>	<i>SG403TX INT</i>	<i>SG603TX INT</i>
<i>Part Name</i>	<i>Part No. SG403TX</i>	<i>Part No. SG603TX</i>
Capacitor	11557	11557
Circuit Breaker - 5A	34331	34331
Exhaust HEPA Filter	38582	38583
Fluorescent Ballast	34160	34160
Fluorescent Lamp	17927	17927
Front Panel, Dress Panel 10 Degree	333752	336572
Ground Fault Interrupter Outlet Indicator Light	34921	34921
Lamp Holder, "Butt-On" Type	34465	34465
Limit Switch	16651	16651
Main Fuss. MDA 20 amp.	31490	31490
Air Flow Monitor (AFM)	38678	38678
Membrane Switch Controller	34699B	34699B
Motor/Blower Controller	304302	304302
Relay, Solid State	1855167	1855167
Sash Balance	38014	
Sash Balance R.H.		38172
Sash Balance L.H.		38173
Side Panel, LH	333X650	333X650
Side Panel, RH	333X651	333X651
Standard Outlet	18231	18231
Supply Blower DD9-4AT	11428	
Supply Blower DD9-6AT		11430
Supply Diffuser	333079	336079
Supply HEPA Filter	37865	37866
Supply Motor 1/4 HP 1625RPM	11619	
Supply Motor 1/3 HP 1625RPM		31676
Switch, Touch-Type	37919	37919
Label, Touchpad, overlay	37998	37998
Transformer	34327	34327
Ultraviolet Lamp	18023	18025
Ultraviolet Lamp Ballast	34160	34160
Ultraviolet Lamp Socket	18088	
Ultraviolet Lamp Socket (plunger type)		20281
Ultraviolet Lamp Socket (fixed type)		20283
Viewscreen Assembly	333066	336066
Viewscreen Frame Assembly	333081	336081