

REPLACEMENT PARTS LIST
STERILGARD III ADVANCE , MODEL SG303 VERTICAL LAMINAR
FLOW CABINET

Part Description	Part #
Capacitor	11557
Circuit Breaker	34331
Exhaust Hepa Filter	39447
Fluorescent Ballast	40093
Fluorescent Lamp	38684
Front Panel	332P605
Ground Fault Interrupter Outlet with Indicator Light	34921
Lamp Holder Base	35072
Lamp Holder Clip	35071
Limit Switch	16651
Mag Gauge	20773
Membrane Switch Controller	34699
Motor/Blower Controller	304227
Relay, Solid State	39636
Sash Balance	39554
Side Panel, LH	332X602
Side Panel, RH	332X601
Standard Outlet	18231
Supply Blower	11417
Supply Diffuser	332P105
Supply Hepa Filter	39861
Supply Motor	39799
Switch, Touch Type	37919
Power Supply	39404
Micro-Controller	39407
Ultraviolet Lamp	18022
Ultraviolet Lamp Ballast	31354
Ultraviolet Lamp Socket	18088
Viewscreen Assembly	332A300
Viewscreen Frame Assembly	332A301