

Replacement Parts List

**STERILGARD® III ADVANCE°, MODEL SG403 & SG603
VERTICAL LAMINAR FLOW CABINET**

Part Name	Part No. SG403	Part No. SG603
Capacitor	11557	11558
Circuit Breaker – 4A & 5A	34331	34337
Exhaust HEPA Filter	37868	37869
Fluorescent Ballast	37868	37869
Fluorescent Light	17927	17927
Front Panel	333068	336155
Ground Fault Interrupter Outlet with Indicator Light	34921	34921
Lamp Holder, “Butt-On” Type	34465	34465
Limit Switch	16651	16651
Magnehelic Gauge	20773	20773
Membrane Switch Controller	34699A	34699A
Motor/Blower Controller	304227	304227
Relay, Solid State	1855167	1855167
Sash Balance	38014	38172, 38173
Side Panel, LH	333152	333152
Side Panel, RH	333153	333153
Standard Outlet	18231	18231
Supply Blower	11429	11432
Supply Diffuser	333079	336079
Supply HEPA Filter	37865	37866
Supply Motor	31676	31679
Switch, Touch-Type	37919	37919
Transformer	34327	34327
Ultraviolet Lamp	18023	18025
Ultraviolet Lamp Ballast	34160	34165
Ultraviolet Lamp Socket	18088	20281, 20283
Viewscreen Assembly	333066	336066
Viewscreen Frame Assembly	333081	336081